



3350 Urbancrest Industrial Drive
 Grove City, Ohio 43123
 Toll Free 800.273.5233
 Fax 614.871.9353
 www.KneeCPM.com

CPM Replacement Parts List 2012

Part Number	Description	Unit Measure
SKIT 1010	Motor Wiring Harness	Each
SKIT 1130	Main Circuit Board (Model 1800)	Each
SKIT 1200	Pendant Cable (Model 1800)	Each
SKIT 1207	Pendant Enclosure Kit (Model 1800)	Each
SKIT 1214	Pendant (Model 1800)	Each
SKIT 1223	Hardware/Fasteners (Model 1800/1850/DC2480)	Each
SKIT 1225	Pendant Start/Stop Switch (Model 1800)	Each
SKIT 1226	Pendant LED Snap Mount (Model 1800)	Each
SKIT 1231	Pendant Velcro Strap (Model 1800) (qty 5)	Pkg
SKIT 1402	Knee Potentiometer Assembly (current—all knee models)	Each
SKIT 1417	Knee Potentiometer Cable Kit	Each
SKIT 1420	Knee Potentiometer Assembly (rectangular)	Each
SKIT 1511	Power Supply, Domestic	Each
SKIT 1600	Foot Assembly Complete (Model 1800/1850/DC2480)	Each
SKIT 1601	Foot Assembly Hardware Kit (upgrade pediatric conversion)	Each
SKIT 1613	Knob Kit (4 male + 2 female knobs)	Each
SKIT 1707 A	Track Seal Kit (qty 4)	Each
Adhesive	Adhesive used with SKIT 1707A Track Seal Kit (1 tube)	Each
SKIT 1727	Delron Spacer & Shoulder Bolt Kit	Each
SKIT 1734	Motor & Gear Box (Model 1800/1850/DC2480)	Each
SKIT 1801	Bed Tube Assembly (qty 1)	Each
SKIT 1804	Cable Velcro Straps (qty 5)	Each
SKIT 1812	Skid Pad Kit (qty 2)	Each
SKIT 1823	Knee Potentiometer Cover Kit w/Hardware (rectangular)	Each
SKIT 1823 A	Knee Potentiometer Covers (rectangular) (qty 3)	Pkg
SKIT 2480-1823	Knee Potentiometer Cover Kit w/Hardware (round)	Pkg.
SKIT 2480-1823 A	Knee Potentiometer Covers (round) (qty 3)	Pkg
SKIT 1825	Range of Motion (ROM) Labels (qty 10)	Pkg
SKIT 1850	Pendant (Model 1850/DC2480)	Each
SKIT 1906	Motor Only (Model 1800/1850/DC2480)	Each
SKIT 2480-1834	Knee CPM Shipping Carton	Each
SKIT 3611	Velcro Kit (foot plate + hip pivot velcro) (qty 6 each)	Each



3350 Urbancrest Industrial Drive
 Grove City, Ohio 43123
 Toll Free 800.273.5233
 Fax 614.871.9353
 www.KneeCPM.com

CPM Replacement Parts List 2012

Part Number	Description	Unit Measure
1223	Pendant Plate Screws (Model 1800)	Each
1228	Pendant Knobs (qty 1)	Each
1301	Thigh Tube Assembly (Right) prior to 2008	Each
1303	Calf Tube Assembly (Right) prior to 2008	Each
1310	Thigh Tube, Pot Assembly (Left) prior to 2008	Each
1312	Calf Tube, Pot Assembly (Left) prior to 2008	Each
1402A	Knee Potentiometer (only) (Model 1800/1850/DC2480)	Each
1507	Main Device On/Off Switch (Model 1800/1850/DC2480)	Each
1603	Foot Tube (Interchangeable left or right)	Each
1604	Pediatric Spacers	Each
1605	Foot Plate (Model 1800/1850/DC2480) DC2480-blue	Each
1607	Carriage Bolt	Each
1608	Foot Plate Screw (small allen head screw inside of foot plate)	Each
1610	White Nylon Spacers	Each
1611	Acorn Nut	Each
1613	Large Female Knob	Each
1614A	Hex Head Bolt 2" (use with female knob on footplate)	Each
1710	Rear End Cap	Each
1711	Front End Cap	Each
1717	Outer Thigh Assembly	Each
1727	Shoulder Bolt	Each
1729	Strut Bearing Spacer	Each
1738	Timing Belt	Each
1740	Motor Coupling Sleeve	Each
1750	Small Male Knob	Each
1827	Top Cover Screw (qty 25)	Each
1850—1850	Pendant Enclosure (Model 1850/DC2480)	Each
1850—2100	PCB Pendant (Model 1850/DC2480)	Each
1850—2102	1850 Main Circuit Board (requires trade-in)	Each
1850—2103	Pendant Overlay (Model 1850/DC2480)	Each
1850—2104	Pendant Cable (Model 1850/DC2480)	Each
2480—1301D	Thigh Tube Assembly (Interchangeable left or right)	Each
2480—1303D	Calf Tube (right side)	Each
2480—1312E	Calf Tube (left/potentiometer side)	Each